

# A new instrument from Metrohm Autolab: VIONIC powered by INTELLO

The only single instrument to offer  $\pm 50V$  compliance voltage,  $\pm 6 A$  current and Electrochemical Impedance Spectroscopy (EIS) up to 10 MHz

Utrecht, March 2021



**VIONIC powered by INTELLO is the new generation instrument created by electrochemical instrumentation innovators Metrohm Autolab. Offering a compact footprint and a sleek and futuristic design never before seen in this category, the VIONIC instrument boasts cutting edge electronics and unique features that improve lab efficiency and safety. The INTELLO software continues Autolab's reputation as creators of powerful software that understands the objective of today's researchers – maximum possibilities, pure discovery.**

## Supplemental features now included at no additional cost

- EIS
- Analog scan
- Second Sense (S2)

VIONIC offers unique features that open up possibilities with your research. **Selectable Floating** making it possible to work with any type of cell setup, independent of the grounding state of the cell components. The **Pure Signal Bridge with Adaptive Cables** expands your connection possibilities while ensuring signal stability.

## Software with time-saving features

VIONIC's INTELLO software accommodates a variety of researchers with workflows for both exploratory and routine measurements. The **Main parameters** feature displays only the parameters you choose; streamlining your workflow and making it easy to perform and delegate routine measurements.

Schedule a demo to see what VIONIC powered by INTELLO can do for you.  
[metrohm.com/demoVIONIC](https://metrohm.com/demoVIONIC)

---

## About Metrohm Autolab

Metrohm Autolab sets the standard for electrochemistry instrumentation. Over 30 years ago, we created the first commercially available digital potentiostat/galvanostat that was completely computer controlled. Today our NOVA software is the most powerful electrochemistry software on the market. With an Autolab potentiostat/galvanostat and NOVA software there are no limits to where your research can go. VIONIC powered by INTELLO our new generation instrument is available from March 1, 2021.

---

## Contact

Aimee Casey, Marketing Communication Manager, Metrohm  
Autolab, +31 655 36 80 92, [aimee.casey@metrohm.com](mailto:aimee.casey@metrohm.com)

