

## 1 Purpose

This document describes the differences between the software version **MagIC Net 2.4 (Build 97)** and the version **MagIC Net 2.4 (97 - 1)**.

## 2 Fixed bugs and problems

### Program part Workplace

During the analysis of solid samples with the Combustion Module's Autosampler (MMS 5000) it was possible that the boat was not correctly lifted.

## 3 Compliance

The current software version does not contain any modifications that affect conformity of **MagIC Net** regarding 21 CFR Part 11 and GAMP.

**Herisau, October 30, 2012**

---



D. Strohm

Executive Vice President, Head of  
R&D



U. Kürsteiner

Head of Quality Management