

1 Purpose

This document describes the differences between the new software version **MagIC Net 2.2 Patch 1** and its predecessor MagIC Net 2.2.

2 Fixed bugs and problems

- Fixed: In rare cases, gradient profiles were not processed correctly.
- Fixed: The column heater of the 881 Compact IC pro heated up to the maximum after a temperature below the current column heater temperature was set.

3 Compliance

The current software version does not contain any modifications regarding 21 CFR Part 11 and GAMP.

Herisau, 16 Feb 2011



D. Strohm

Vice President, Head of R&D