

Mira Cal Software

Release Notes

8.0105.8008EN / 2018-03-28



Metrohm AG
Ionenstrasse
CH-9100 Herisau
Switzerland
Phone +41 71 353 85 85
Fax +41 71 353 89 01
info@metrohm.com
www.metrohm.com

Mira Cal Software

3.1.20

Release Notes

8.0105.8008EN /
2018-03-28

Technical Communication
Metrohm AG
CH-9100 Herisau
techcom@metrohm.com

This documentation is protected by copyright. All rights reserved.

This documentation has been prepared with great care. However, errors can never be entirely ruled out. Please send comments regarding possible errors to the address above.

1 Purpose

This document describes the fixed bugs and problems as well as improvements in the following new software versions in comparison with the previous versions.

- New: 3.1.20
- Old: 3.1.18



2 New features

None



3 Improvements

None



4 Fixed bugs and problems

Topic	Description
Firmware updates	Fixed an issue with updating firmware. The previous version reported: "Firmware Update Error... The directory is not empty"
Batch and container fix	There was an error when there was no Lot ID present, other fields would not populate.
Sign level 2 window title	The sign level 2 window title incorrectly read: "Sign Level 1"
Spelling	Capitalized "ID" everywhere.



5 Compliance

The new software version fixes an issue with data integrity and thus improves MiraCal compliance with 21 CFR Part 11 and GMP.