

Mira DS firmware package 8.0.0.98

Release Notes

8.0926.8002EN / 2018-08-24





Metrohm AG
Ionenstrasse
CH-9100 Herisau
Switzerland
Phone +41 71 353 85 85
Fax +41 71 353 89 01
info@metrohm.com
www.metrohm.com

Mira DS firmware package 8.0.0.98

8.0.0.98

Release Notes

8.0926.8002EN /
2018-08-24

Technical Communication
Metrohm AG
CH-9100 Herisau
techcom@metrohm.com

This documentation is protected by copyright. All rights reserved.

This documentation has been prepared with great care. However, errors can never be entirely ruled out. Please send comments regarding possible errors to the address above.

1 Purpose

This document describes the fixed bugs and problems as well as improvements in the following new firmware version in comparison with the previous version.

- New version: 8.0.0.98
- Old version: 8.0.0.91



2 New features

| Topic | Description |
|--------------------|--|
| New hardware | The "ID Kit" and SERS test kit were added. |
| Library management | Library management can be enabled and disabled. |
| Reprocess samples | Reprocess samples by using operating procedure parameters against enabled libraries. |
| Added languages | Czech, Hungarian and Turkish were added. |

