



# Compatibility of *tiamo*<sup>TM</sup> versions with Metrohm devices

Device	<i>tiamo</i> <sup>TM</sup> 1.0	<i>tiamo</i> <sup>TM</sup> 1.1	<i>tiamo</i> <sup>TM</sup> 1.2	<i>tiamo</i> <sup>TM</sup> 1.2.1	<i>tiamo</i> <sup>TM</sup> 1.3	<i>tiamo</i> <sup>TM</sup> 2.0	<i>tiamo</i> <sup>TM</sup> 2.1	<i>tiamo</i> <sup>TM</sup> 2.2	<i>tiamo</i> <sup>TM</sup> 2.3	<i>tiamo</i> <sup>TM</sup> 2.4	<i>tiamo</i> <sup>TM</sup> 2.5	<i>tiamo</i> <sup>TM</sup> 3.0
<b>Titrimos</b>	702 SM Titrimo	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	716 DMS Titrimo	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	718 STAT Titrimo	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	719 S Titrimo	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	720 KFS Titrimo	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	721 Net Titrimo	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	736 GP Titrimo	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	751 GPD Titrimo	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	758 KFD Titrimo	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	784 KFP Titrimo	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	785 DMP Titrimo	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	794 Basic Titrimo	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	795 KFT Titrimo	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	798 MPT Titrimo	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	799 GPT Titrimo	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
<b>Titrandos</b>	808	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	809	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	835	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	836	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	841	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	842	-	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	851	-	-	-	-	-	-	✓	✓	✓	✓	✓
	852	-	-	-	-	-	-	✓	✓	✓	✓	✓
	855	-	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	857	-	-	✓	✓	✓	✓	✓	✓	✓	✓	✓
	888	-	-	-	-	✓	✓	✓	✓	✓	✓	✓
	890	-	-	-	-	✓	✓	✓	✓	✓	✓	✓
	901	-	-	-	-	-	✓	✓	✓	✓	✓	✓
	902	-	-	-	-	-	✓	✓	✓	✓	✓	✓
	904	-	-	-	-	-	✓	✓	✓	✓	✓	✓
	905	-	-	-	-	-	✓	✓	✓	✓	✓	✓
	906	-	-	-	-	-	✓	✓	✓	✓	✓	✓
907	-	-	-	-	-	✓	✓	✓	✓	✓	✓	



Device	<i>tiamo</i> <sup>™</sup> 1.0	<i>tiamo</i> <sup>™</sup> 1.1	<i>tiamo</i> <sup>™</sup> 1.2	<i>tiamo</i> <sup>™</sup> 1.2.1	<i>tiamo</i> <sup>™</sup> 1.3	<i>tiamo</i> <sup>™</sup> 2.0	<i>tiamo</i> <sup>™</sup> 2.1	<i>tiamo</i> <sup>™</sup> 2.2	<i>tiamo</i> <sup>™</sup> 2.3	<i>tiamo</i> <sup>™</sup> 2.4	<i>tiamo</i> <sup>™</sup> 2.5	<i>tiamo</i> <sup>™</sup> 3.0
<b>Coulometers</b>	756 KF Coulometer	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	831 KF Coulometer	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
<b>Sample Changers</b>	730 Sample changer	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	774 Oven Sample Changer	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	778 Sample Processor	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	789 Robotic SP XL	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	814 USB Sample Processor	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	815 Robotic USB Sample Processor	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	855 Robotic Titrosampler	–	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	864 Robotic Balance Sample Processor	–	–	–	–	–	✓	✓	✓	✓	✓	✓
874 Oven Sample Processor	–	–	–	–	–	✓	✓	✓	✓	✓	✓	
<b>Meters</b>	712 Conductometer	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	856 Conductivity Module	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
<b>Miscellaneous</b>	846 Dosing Interface	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	854 iConnect	–	–	✓	✓	✓	✓	✓	✓	✓	✓	✓
	859 Titrotherm	–	–	–	–	–	–	✓	✓	✓	✓	✓
	867 pH module	–	–	–	–	–	✓	✓	✓	✓	✓	✓
	Avantes spectrometer	–	–	–	–	–	✓	✓	✓	✓	✓	✓
	Barcode readers	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓