

## 1 Declaration of Conformity

### *Description*

Database software for the storage, visualisation and recalculation of titration data. Depending on the analysis device data from PC/LIMS-reports are imported via USB-stick, CompactFlash card, USB lab link or a RS-232 interface and transformed into a database readable form. Client/server functions are possible. Software runs under Windows™ 2000, Windows™ XP, Windows™ Vista or Windows™ 7. Dialog in German, English.

**6.6063.112****tiBase 1.1 full**

Local server version with 1 license

**6.6063.113****tiBase 1.1 multi**

Client/server version with 3 licenses

### *Quality system*

The software mentioned above was developed in accordance with the requirements of the ISO 9001:2008 quality system regarding the design, testing and servicing of Metrodata software.

### *Technical specifications*

The technical specifications are documented in the installation manual.

### *Software validation*

The software was validated with respect to functionality, analytical performance and accuracy of results. The software functions are documented in the online help.

**Herisau, 29 June 2012**



D. Strohm

Executive Vice President, Head of  
R&D



U. Kürsteiner

Head of Quality Management