



## Compatibility of *Vision Air* versions with Windows™ operating systems

	Windows XP (SP1, SP2, SP3)	Windows Vista	Windows 7	Windows 8	Windows 10
<i>Vision Air</i> 1.1	-	-	√	√	-
<i>Vision Air</i> 2.0	-	-	√	-	√

only 64-bit operating systems are supported

	Windows Server 2008 R2 32&64 bit	Windows Server 2012 64 bit	Windows Server 2012 R2 64 bit
<i>Vision Air Server</i> 1.1	√	-	-
<i>Vision Air Server</i> 2.0	-	√	√