



OMNIS meets CUBIS® II

Seamless integration
for enhanced
efficiency, compliance,
and data integrity

PEOPLE
YOU
CAN
TRUST

Cubis® II laboratory balance is part of OMNIS



DISCOVER THE BENEFITS FROM THE INTEGRATION OF CUBIS® II INTO THE OMNIS WORLD

	Cubis® II	Non-Cubis® Balances
Transfer of Weighing Values	Yes	Yes
Transfer of Metadata and Audit Trail	Yes	Limited
Plug-and-Play Intergration	Yes	No
Balance Sharing with Multiple OMNIS Systems	Yes	No
Seamless Workflow Support	Yes	Limited
Operation Directly on the Balance's Display	Yes	No
Centralized User Management	Yes	No

IN TODAY'S LABORATORIES, PRECISION, EFFICIENCY, AND REGULATORY COMPLIANCE ARE CRITICAL TO SUCCESS.

The powerful combination of OMNIS by Metrohm and the Sartorius Cubis® II premium laboratory balance delivers an integrated solution for automated sample preparation – from accurate weighing to titration.

This seamless workflow eliminates manual data transfer, saves valuable time, and supports full compliance with 21 CFR Part 11 and EU Annex 11 requirements for data integrity. Thanks to its modular system architecture, the solution adapts to your lab's evolving needs.

EFFORTLESS, SECURE DATA HANDLING

More than numbers: Cubis® II sends encrypted results and metadata to OMNIS, ensuring data integrity and compliance.

ONE SYSTEM, FULL CONTROL

OMNIS centrally manages user access and balance activity with a real-time Audit Trail – fully compliant with FDA 21 CFR Part 11 and Annex 11, effortlessly.

SMART WEIGHING AT THE BALANCE

Start sample series from Cubis® II without opening OMNIS. Visual guides support weighing, while multiple OMNIS clients share one balance for efficient workflows.

