



Pure Electrocatalysis

VIONIC powered by INTELLO offers the highest combined specifications in one single instrument.

HIGHLIGHTS

- **Pure Efficiency:** VIONIC and the INTELLO software are bursting with time-saving features that optimize any workflow.
- **Pure Versatility:** With unmatched standard specifications, VIONIC is the instrument of choice for your electrochemical applications.
- **Pure Safety:** Smart hardware and software safety features protect your cell, your lab, and your data.
- **Pure Discovery:** VIONIC powered by INTELLO offers a unique combination of features that observes electrochemical processes, in real time with no gaps or missed information: complete data, pure discovery.

VIONIC powered by INTELLO offers the highest combined specifications in one single instrument

Maximum measured current	± 6 A, over 11 current ranges
Dual* mode compliance voltage	± 50V
Maximum** applied potential	± 10V
Second sense	Yes, additional voltage measurement in the cell
EIS frequency range	from 10 μ Hz to 10 MHz
True analog scan	from 50 μ V/s to 100 V/s
Lowest sampling interval (acquisition) time (i, E, S2)	1 μ s

*High accuracy or high voltage **Maximum compliance mode

Dedicated to research: metrohm.com/electrochemistry



OTHER FEATURES

- **Seamless** measurements in real-time, with no gaps or missed information providing the most complete data.
- **Untethering** gives you the freedom to use your computer for other activities while your experiment runs.
- Up to **10 million data points** can be stored, more than enough for typical electrochemical measurements.

METROHM AUTOLAB

The research and development of VIONIC powered by INTELLO was based upon over 30 years of customer insight. Each feature was created to meet the requirements of electrochemical research and improve your day-to-day discoveries.

Metrohm Autolab is an ISO 9001 certified company.



Metrohm Autolab provides an industry-leading 3 year warranty on all instruments and accessories.