

Unitrode with plug-in head U



NEW

The proven Unitrode with Pt 1000 temperature sensor now also comes with a removable electrode cable



The advantages at a glance

- **Cost advantage:** Lower replacement cost than for an electrode with fixed cable.
- **Simple storage:** The electrode is easy to store between measurement series. There is no cable in the way, as it can be stored separately or remains on the instrument.
- **Faster electrode changing – easy to assemble even in automated systems:** The cable remains permanently installed in the titration system; only the sensor on the titration head has to be changed.
- **Metrohm's contribution to environmental protection:** Reduction in cable waste.
- **Reliability:** Thanks to gold-plated contacts, the unique multi-pin connector for lab electrodes offers a high level of signal reliability.

Ordering information

6.0258.600	LL Unitrode Pt 1000 (plug-in head U)
6.2104.600	Electrode cable with plug-in head U / connector F + 2 x B, 1 m
6.2104.610	Electrode cable with plug-in head U / connector F + 2 x B, 2 m

Visit us at www.metrohm.com