

# Reliability for Routine pH Measurements



## NEW

The proven Ecotrode now available as a maintenance-free gel electrode with aging indicator



Still within the permissible operating range: The decolorized zone of the gel gradually spreads up to the mark with use.

## Highlights

- Ideal for **routine measurements** with similar types of samples
- **Maintenance-free gel electrolytes:** Refilling or checking the electrolyte is not necessary
- **Reliability through aging indicator:** The white gel loses its color with use.  
As long as the decolorized zone is below the marking, accurate results are guaranteed.

### Ordering information

|            |                                                               |
|------------|---------------------------------------------------------------|
| 6.0221.100 | Ecotrode Gel (plug-in head G)                                 |
| 6.0221.600 | Ecotrode Gel (plug-in head U)                                 |
| 6.2104.020 | Electrode cable with plug-in head G / connector F             |
| 6.2104.600 | Electrode cable with plug-in head U / connector F, 2 x B, 1 m |

Visit us at [www.metrohm.com](http://www.metrohm.com)