

Unitrode easyClean



High-performance pH sensor with diaphragm cleaning at the press of a button

NEW: The proven all-rounder with Pt1000 temperature sensor now also with easyClean diaphragm – truly a clean affair!




The benefits at a glance

- For universal use, even in difficult samples
- High temperature resistance and very low alkali error
- Simple, contact-free cleaning of the diaphragm
- Reduced immersion depth through optimized membrane shape
- Including Pt1000 temperature sensor
- Optional: Idrolyte as electrolyte for temperatures of 80–100 °C
- For Titrino plus, Ti-Touch and Titrande family instruments

Ordering information

6.0260.010	Unitrode easyClean, 1.2 m fixed cable, plug F + 1 x 2 mm
6.0260.020	Unitrode easyClean, 2 m fixed cable, plug F + 1 x 2 mm
6.2308.040	Idrolyte, 250 mL
6.2308.020	Electrolyte 3 mol/L KCl (250 mL)
6.2323.000	Storage solution
6.2325.000	pHit kit



www.metrohm.com

 **Metrohm**