



900 Touch Control

Titration was never easier

PEOPLE
YOU
CAN
TRUST

 **Metrohm**

Titration at a single touch

Powerful as a PC, smaller than a netbook and a brilliant touchscreen – that is the 900 Touch Control. This benchtop control module makes titration easier, safer and faster than ever before. These are the 900 Touch Control's major functions and user benefits:

USB INTERFACE

USB interface for printer, keyboard or a mouse. Simply plug in a USB stick to save your methods and determinations.

TOUCH AND TITRATE

Up to 14 methods can each be linked to a favorite icon on the touchscreen. You start your standard methods at a single touch. Titration cannot be easier!



FDA COMPLIANCE

The 900 Touch Control complies with FDA-regulation 21 CFR, Part 11.

ETHERNET & LIMS

Ethernet interface for direct access to your intranet. Print reports on a network printer or save them on a network storage. Send your results directly to your LIMS or store them in the 900 Touch Control's archive.

MULTILINGUAL USER DIALOG

Besides German, English, French and Spanish, the 900 Touch Control communicates also in Chinese, Russian, and Japanese.

AUTOMATIC ELECTRODE TEST

The 900 Touch Control comes with an automatic electrode test. Damaged electrodes are identified immediately; you can always rely on correct results.

AUTOSTART FUNCTION

In Karl Fischer titration it can happen all too fast that sample is injected into the titration cell without the start button for the titration having been pressed first. In this case the sample is conditioned away and the determination has to be repeated. This cannot happen with the 900 Touch Control as the device starts the titration automatically once sample is added.



PDF GENERATOR

Use the 900 Touch Control's PDF generator to create forge-proof reports – without a PC!

ORDERING INFORMATION

Instrument

2.900.0010 900 Touch Control

The 900 Touch Control is compatible with:

808, 809, 835, 836, 841, 842, 851, 852, 857, 888, 890, 901, 902, 904, 905, 906, 907 Titrandos

846 Dosing Interface

856 Conductivity Module

867 pH Module

814 USB Sample Processor

815 Robotic USB Sample Processor

815 Robotic Soliprep

855 Robotic Titrosampler

855 Robotic Analyzers

885 Compact Oven Sample Processor
