



Compatibility of *Vision* versions with Windows™ operating systems

	Windows 2000	Windows XP (SP1/SP2/SP3)	Windows Vista	Windows 7	Windows 8	Windows 10	Windows 11
Vision 4.0.1–4.0.2	√	√	√	√	√*	–	–
Vision Multi 4.0.1–4.0.2	√	√	√	√	√*	–	–
Vision Pharma 4.0.1–4.0.2	√	√	√	√	√*	–	–

* only 64-bit operating systems are supported

Vision 4.0.3	√	√	√	√	√	–	–
Vision Multi 4.0.3	√	√	√	√	√	–	–
Vision Pharma 4.0.3	√	√	√	√	√	–	–

Vision 4.1–4.1.1.156	–	–	√	√	√	√	–
Vision Multi 4.1–4.1.1.156	–	–	√	√	√	√	–
Vision Pharma 4.1–4.1.1.156	–	–	√	√	√	√	–
Vision Multi Pharma 4.1–4.1.1.156	–	–	√	√	√	√	–

Vision 4.1.1.189–4.1.1.191	–	–	–	√	√	√	–
Vision Multi 4.1.1.189–4.1.1.191	–	–	–	√	√	√	–
Vision Pharma 4.1.1.189–4.1.1.191	–	–	–	√	√	√	–
Vision Multi Pharma 4.1.1.189–4.1.1.191	–	–	–	√	√	√	–

Vision 4.1.1.238 or higher	–	–	–	√	√	√	√
Vision Multi 4.1.1.238 or higher	–	–	–	√	√	√	√
Vision Pharma 4.1.1.238 or higher	–	–	–	√	√	√	√
Vision Multi Pharma 4.1.1.238 or higher	–	–	–	√	√	√	√