

## **VA Database 2.1**

**Beilage 8.746.8037  
zur Gebrauchsanweisung 8.746.8003**

### **Neue Programmversion**

Das Programm «VA Database 2.1» kann sowohl Daten der neuen Programmversion 5.746.0101 des VA Processors 746 mit geändertem Datumsformat wie auch Daten der alten Version 5.746.0100 verarbeiten. Ansonsten ist es identisch mit der Version «VA Database 2.0».

### **Bestellbezeichnungen**

**VA Database 2.1 (Vollversion) ..... 6.6016.113**  
enthält Software und Dongle 6.2145.030

**VA Database 2.1 (Demoversion/Update)..... 6.6016.118**  
enthält Software ohne Dongle 6.2145.030

**Dongle ..... 6.2145.030**

## **VA Database 2.1**

**8.746.8037 Annex  
to the 8.746.8003 Instructions for Use**

### **New program version**

The program «VA Database 2.1» can process data of the new program version 5.746.0101 of the 746 VA Processor with changed format of date as well as data of the old version 5.746.0100. Otherwise it is identical to the version «VA Database 2.0».

### **Ordering designations**

**VA Database 2.1 (Full version)..... 6.6016.113**  
including Software and 6.2145.030 Dongle

**VA Database 2.1 (Demo version/Update)..... 6.6016.118**  
including Software without 6.2145.030 Dongle

**Dongle..... 6.2145.030**