

VA Database 2.2

**Beilage 8.746.8047
zur Gebrauchsanweisung 8.746.8003**

Neue Programmversion

Das Programm «VA Database 2.2» kann sowohl Daten der neuen Programmversion 5.746.0101 des VA Processors 746 mit geändertem Datumsformat wie auch Daten der alten Version 5.746.0100 verarbeiten. Neu besteht die Möglichkeit, Daten automatisch in eine txt-Datei zu exportieren (Details siehe Help). Ansonsten ist es identisch mit der Version «VA Database 2.0».

Bestellbezeichnungen

VA Database 2.2 (Vollversion)	6.6016.123
enthält Software und Dongle 6.2145.030	
VA Database 2.2 (Demoversion/Update)	6.6016.128
enthält Software ohne Dongle 6.2145.030	
Dongle	6.2145.030

VA Database 2.2

**8.746.8047 Annex
to the 8.746.8003 Instructions for Use**

New program version

The program «VA Database 2.2» can process data of the new program version 5.746.0101 of the 746 VA Processor with changed format of date as well as data of the old version 5.746.0100. A new feature is the possibility to export automatically data into a txt file (details see Help). Otherwise it is identical to the version «VA Database 2.0».

Ordering designations

VA Database 2.2 (Full version)	6.6016.123
including Software and 6.2145.030 Dongle	
VA Database 2.2 (Demo version/Update)	6.6016.128
including Software without 6.2145.030 Dongle	
Dongle	6.2145.030