

1 Corrigendum (Online-Hilfe 8.101.8025DE)

Kapitel 1.2 Anwendung

Alt tiBase ist ein Datenbankprogramm zur Verwaltung von Daten auf der Basis von PC/LIMS-Reports. Die Reports von den entsprechenden Metrohm-Geräten (Titrino, Titrino Plus, Titrande und Dosimat plus) können in tiBase eingelesen werden.

Neu tiBase ist ein Datenbankprogramm zur Verwaltung von Daten auf der Basis von PC/LIMS-Reports. Die Reports von Titrino plus, Compact Titrosampller, PC Control 6.0 und Touch Control 5.840.0150 können in tiBase eingelesen werden.

Kapitel 1.4 Geräteeinbindung

Alt Von folgenden Geräten können PC/LIMS-Reports in tiBase eingelesen werden:

- Titrando
808, 809, 835, 836, 841, 842, 855, 857, 888, 890, 901, 902, 904, 905, 906, 907
- Titrino
702, 716, 718, 719, 720, 721, 736, 751, 758, 784, 785, 794, 795, 798, 799
- Titrino plus
848, 862, 870, 877
- Dosimat plus
865, 876

Neu PC/LIMS-Reports können in tiBase eingelesen werden von:

- Titrino plus
848, 870, 877
- 862 Compact Titrosampller
- Titrando
808, 809, 835, 836, 841, 842, 855, 857, 888, 890, 901, 902, 904, 905, 906, 907
Steuerung mit Touch Control 5.840.0150 oder PC Control 6.0

2 Corrigendum (online help 8.101.8025EN)

Chapter 1.2 Application

Old tiBase is a database program for administrating data on the basis of PC/LIMS reports. The reports from the corresponding Metrohm instruments (Titrino, Titrino Plus, Titrande and Dosimat plus) can be read into tiBase.

New tiBase is a database program for administrating data on the basis of PC/LIMS reports. The reports from Titrino plus, Compact Titrosampler, PC Control 6.0 and Touch Control 5.840.0150 can be read into tiBase.

Chapter 1.4 Connecting devices

Old PC/LIMS reports of the following instruments can be read into tiBase:

- Titrande
808, 809, 835, 836, 841, 842, 855, 857, 888, 890, 901, 902, 904, 905, 906, 907
- Titrino
702, 716, 718, 719, 720, 721, 736, 751, 758, 784, 785, 794, 795, 798, 799
- Titrino plus
848, 862, 870, 877
- Dosimat plus
865, 876

New PC/LIMS reports can be read into tiBase from:

- Titrino plus
848, 870, 877
- 862 Compact Titrosampler
- Titrande
808, 809, 835, 836, 841, 842, 855, 857, 888, 890, 901, 902, 904, 905, 906, 907
Control by Touch Control 5.840.0150 or PC Control 6.0