

## 1 Purpose

This document describes the differences between the new software version **Metrohm IC Driver 2.0 for Empower™** and its predecessor **Metrohm IC Driver 1.1 for Empower™**.

This release of **Metrohm IC Driver 2.0 for Empower™** replaces all previous releases.

## 2 Improvements

1. Integration of device family **940 Professional IC Vario**.  
(version: 2.940.2400, 2.940.2500, 2.940.1440, 2.940.1540)
2. Integration of device family **942 Extension Module Vario**.  
(version: 2.942.0040, 2.942.1060)
3. Integration of device family **889 IC Sample Center**.  
(version: 2.889.0010, 2.889.0020)
4. Integration of device **945 Professional Detector Vario – Amperometry**.  
(version: 2.945.0020)
5. Integration of device **Remote Box MSB**.  
(version: 6.2148.010)

## 3 Fixed bugs and problems

none

## 4 Known problems for resolution in next release


none

Herisau, October 13, 2015

---



P. Hunziker  
Vice President,  
Head of Development



M. Matter  
Head of Quality Management