

1 Purpose

This document describes the fixed bugs as well as improvements in the following new software version in comparison with the previous version.

- New version: 4.1.11
- Old version: 4.1.1

2 New features

None

3 Improvements

Topic	Description
Added pass/fail to audit trail	Added pass/fail result to the details of the sample audit trail.
Audit trail entries for shutdown	Added audit trail coverage for shutdown events.

4 Fixed bugs

Topic	Description
Result consistency	Fixed result processing by truncating values for consistency between instrument and software.
Filename adjustment to prevent overwrite on export	Added fractional seconds to the filenames to prevent identical filenames from occurring.

5 Known issues

None.

6 Compliance

The current software version does not contain any modifications that affect conformity of Mira Cal P regarding 21 CFR Part 11 and GAMP.