

1 Purpose

This document describes the fixed bugs in the following new software version in comparison with the previous version.

- New version: 4.2.200
- Old version: 4.2.100

2 New features

None.

3 Improvements

None.

4 Fixed bugs

Topic	Description
SST sync fails in other language	System Suitability Tests no longer causes a synchronization error when Mira Cal P was set to a language other than English.
Database conversion > 4 GB	Database conversion process is fixed so that databases > 4 GB will not crash. <i>Applies to conversion from 4.1.X to 4.2.200</i>
LastSST field is populated in XML	The LastSST field is correctly populated now when saving sample data as XML or when building Custom Reports.

5 Known issues

None.

6 Compliance

The current firmware version does not contain any modifications that affect conformity of Mira Cal regarding 21 CFR Part 11 and GAMP.