

## 1 Purpose

This document describes the fixed bugs in the following new software version in comparison with the previous version.

- New version: 4.2.200
- Old version: 4.2.100

## 2 New features

Topic	Description
Windows 11	Compatible with Windows 11.  System requirements (in addition to the Windows 11 system requirements): <ul style="list-style-type: none"><li>▪ Storage: 100 GB minimum, 500 GB recommended</li><li>▪ USB connectors: USB 3.0 Note: Only use the provided Metrohm USB cable (6.021.08010), do not use third party USB cables. The usage of a powered USB Hub is recommended.</li><li>▪ Mouse and keyboard: Required, USB, PS/2 style, or wireless</li><li>▪ Screen resolution 1024 × 768 or higher</li><li>▪ System backup: Network or local backup required for data archiving.</li></ul>

## 3 Improvements

None.

## 4 Fixed bugs

Topic	Description
SST sync fails in other language	System Suitability Tests no longer causes a synchronization error when Mira Cal P was set to a language other than English.

Topic	Description
Database conversion > 4 GB	Database conversion process is fixed so that databases > 4 GB will not crash.  <i>Applies to conversion from 4.1.X to 4.2.200</i>
<b>LastSST</b> field is populated in XML	The <b>LastSST</b> field is correctly populated now when saving sample data as XML or when building Custom Reports.

**5 Known issues**

None.

**6 Compliance**

The current firmware version does not contain any modifications that affect conformity of Mira Cal regarding 21 CFR Part 11 and GAMP.