

1 Purpose

This document describes the fixed bugs and problems as well as improvements in the following new software version in comparison with the previous version.

- New version: 1.2.55
- Old version: 1.2.50

2 New features

Topic	Description
New language packs	id – Indonesian ja – Japanese ko – Korean pl – Polish ru – Russian th – Thai vi – Vietnamese

3 Improvements

Topic	Description
Libraries	Supports the following library: <ul style="list-style-type: none">▪ 66071634 MIRA Illicit & General Chemicals Library

4 Fixed bugs

Topic	Description
Chinese translations	Fixed translations for various topics.
German translations	Fixed translations for various topics.
cs, de, es, fr, hu, it, pt, tr, zh language packs	If the service recovery tool was used, the language packs caused an issue that forced the device to reboot after every 5 time it was powered up.

5 Known issues

- No support for Arabic language.
- Some language packs have errors or missing strings in the UI.
- One cannot view XTR processed data and baselined data for the same spectrum at the same time. One cannot view baselined data of an XTR processed spectrum.