

## 1 Purpose

This document describes new functions, the changes, improvements and fixed bugs and problems in the following firmware version:

- 57.1007.0011

in comparison with the firmware version:

- 57.1007.0009

## 2 New functions

The following new functions have been integrated:

- Changing the password for the **Expert** dialog type
- Factory reset
- User interface: 10 new languages added:
  - French
  - Italian
  - Japanese
  - Korean
  - Polish
  - Portuguese
  - Russian
  - Spanish
  - Thai
  - Vietnamese

## 3 Improvements

The following improvements were made:

- The Eco Dosimat user interface is aligned with the Eco Titrator user interface concept.
- Firmware optimizations for the final instrument check in the production.

