

1 Purpose

This document describes new functions, the changes, improvements and fixed bugs in the following firmware version:

- 57.1007.0017 for OS version 007

In comparison with the firmware version:

- 57.1007.0016 for OS version 007

2 New functions

The following new functions have been integrated:

- Initial adjustment for servo motor at time of device production
- Detection of main board version in relation to initial adjustment

3 Improvements

The following improvements were made:

- Optimizations for tests in production

4 Fixed bugs

The following bugs were fixed:

- Sporadic printing errors (TITR-22037)

5 Known issues

The following bugs are known:

- The PDF report also displays parameters that have not been activated.