

## 1 Purpose

This document describes new functions, changes, and improvements in the following firmware version:

- 57.1007.0019 for OS version 008

In comparison with the firmware version:

- 57.1007.0018 for OS version 008

## 2 New functions

The following new functions have been integrated:

- Various new functions for service and production.

## 3 Improvements

The following improvements were made:

- None

## 4 Fixed bugs

The following bugs were fixed:

- None

## 5 Known issues

The following bugs are known:

- If an attempt is made to output the parameters as a PDF or via a Post-script printer directly after a factory reset, then the report will be empty.