

## 1 Purpose

This document describes new functions, the changes, improvements and fixed bugs in the following firmware version:

- 57.1008.0015 for OS version 006
- in comparison with the firmware version:
- 57.1008.0014 for OS version 006

## 2 New functions

The following new functions have been integrated:

- Export of CSV files (TITR-18026).

## 3 Improvements

The following improvements were made:

- Occurrence of **Piston blocked** error message was regulated.

## 4 Fixed bugs

The following bugs were fixed:

- Temporary freezing after pre-dosing (TITR-18647).
- Erroneous calculations (TITR-18812)
- Titration in SET mode stops before the endpoint is reached (TITR-19339).

## 5 Known bugs

The following bugs are known:

- The PDF report also displays parameters that have not been activated (TITR-16729).
- Sporadic printing errors (TITR-17204)