

1 Purpose

This document describes the new features, the improvements as well as the fixed bugs and problems in the following new software:

- Version 1.1

in comparison with the software:

- Version 1.0

2 New features

Topic	Description
Service area	The software now features a service area for the diagnosis and error detection in the measuring instrument.
Library	The update of S.T. Japan libraries was implemented.

3 Improvements

Topic	Description
Numerical values	The display of numerical values with their decimal places has been standardized.
Data transfer / data evaluation	The data transfer from the standalone application of the instrument has been improved.
Library	The actualization of custom libraries has been improved. The data format of the S.T. Japan libraries was expanded.
Evaluations	Reports have been improved, so that they can now be used for several samples.
User interface	Several details on the user interface have been improved.

4 Fixed bugs and problems

Topic	Description
Report	The report error when making evaluations with no matches in the library was eliminated. The error when giving the report file an already existing name was eliminated.
Time synchronization	The faulty time synchronization between instrument and software in different time formats was corrected.


5 Known errors

Topic	Description
none	none

Herisau, May 28, 2015



P. Hunziker
Vice President,
Head of Development



M. Matter
Head of Quality Management