

## 1 Purpose

This document describes new features and fixed bugs in the following new software version in comparison with the previous version.

- New version: 4.1.1.54
- Old version: 4.1.1.51

## 2 New features

Topic	Description
New Standard Set	Standard Set with order number 6.7450.130 can be used for instrument calibration, wavelength certification, and photometric tests.

---

## 3 Fixed Bugs

Topic	Description
Vision Air File Import	<p>In certain cases constituent values were incorrectly imported. This has been fixed.</p> <p>Vision crashed when importing an XML file from Vision Air, which contained a timeout entry from Smart Probe. This has been fixed.</p>
Vision C++ Redistributable	In some cases the 32 bit version of Vision C++ Redistributable would not be installed on 64 bit systems, which caused the Vision installation to fail. This has been fixed.

---

## 4 Compliance

The current software version does not contain any modifications that affect conformity of Vision regarding 21 CFR Part 11 and GAMP.