

1 Purpose

This document describes the improvements in the following new firmware version in comparison with the previous versions:

- 5.888.0023


2 Improvements

Topic	Description
KFT Ipol, KFT Upol	The 888 Titrande now supports the titration modes KFT Ipol and KFT Upol (in conjunction with tiamo 2.5).


3 Compliance

The current firmware version does not contain any modifications that affect the conformity of the instrument regarding 21 CFR Part 11.

Herisau, September 15, 2015



P. Hunziker
Vice President,
Head of Development



M. Matter
Head of Quality Management