

1 Purpose

This document describes the new features, the improvements as well as the fixed bugs and problems in the following new firmware versions:

- 5.912.0013
- 5.913.0013
- 5.914.1013
- 5.914.2013

in comparison with the firmware versions:

- 5.912.0012
- 5.913.0012
- 5.914.1012
- 5.914.2012

2 New features

Topic	Description
none	none

3 Improvements

Topic	Description
Energy supply	Charging of (almost) flat or low batteries was improved. The battery status is now indicated more accurately.

4 Fixed bugs and problems

Topic	Description
none	none

5 Known bugs

Topic	Description
none	none

Herisau, June 15, 2017



P. Hunziker
Vice President,
Head of Development



M. Saba
Head of Quality Management