



OMNIS Reagent Organizer Liquid Control

2.1056.1200

The OMNIS Reagent Organizer Liquid Control serves as Reagent Organizer and provides the illumination of the interior in the OMNIS IC System. In addition, it offers a software-integrated leak detection.

The OMNIS Reagent Organizer Liquid Control is always located on the topmost OMNIS Extension Module and also forms the lid of the OMNIS IC System. A leak sensor is installed at the drain of the OMNIS Reagent Organizer Liquid Control. The lighting integrated into the OMNIS Reagent Organizer Liquid Control illuminates the interior of the OMNIS IC System as soon as at least one door is opened.

The OMNIS Reagent Organizer Liquid Control offers space for up to 6 bottles and can be equipped with additional accessories, as for example bottle holder, bottle fixation, volume sensor, magnetic stirrer or waste container with external canister sensor.

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 2.1056.1200

Order Number	Description	Quantity
1.1056.1200	OMNIS Reagent Organizer Liquid Control The OMNIS Reagent Organizer Liquid Control serves as Reagent Organizer and provides the illumination of the interior in the OMNIS IC System. In addition, it offers a software-integrated leak detection. The OMNIS Reagent Organizer Liquid Control is always located on the topmost OMNIS Extension Module and also forms the lid of the OMNIS IC System. A leak sensor is installed at the drain of the OMNIS Reagent Organizer Liquid Control. The lighting integrated into the OMNIS Reagent Organizer Liquid Control illuminates the interior of the OMNIS IC System as soon as at least one door is opened. The OMNIS Reagent Organizer Liquid Control offers space for up to 6 bottles and can be equipped with additional accessories, as for example bottle holder, bottle fixation, volume sensor, magnetic stirrer or waste container with external canister sensor.	1 PCS
6.02102.010	Cable MDL PL/SO 0.5 m Connecting cable for OMNIS instruments with MDL connector (0.5 m)	1 PCS

Optional accessories

6.01602.100

Bottle cap GL 45, 2 x M6, UNF 10/32

Bottle cap GL 45 for eluent, regenerant, reagent, and rinsing bottles, with connectors for the adsorber tube and the aspiration tubing.



6.01602.110

Bottle cap S 40, 2x M6, UNF 10 /32

Cap S 40 for 10 L containers.



6.01608.000

Bottle PP, 2 L, GL 45

Bottle 2 L (PP) designed for OMNIS Reagent Organizer, with volume indication for eluent, regenerant, reagent, and rinsing solutions.



6.02009.000

**Volume sensor holder,
container 10 L**

Holder for 10 L container with integrated volume sensor.



6.02050.000

Bottle holder, RO

Bottle holder for Eluent bottle PP 2 L on OMNIS IC Reagent Organizer.



6.02050.010

Bottle fixation RO

Bottle fixation for Eluent bottle PP 2 L on OMNIS IC Reagent Organizer.



6.02050.100

Bottle holder RO, universal

Bottle holder for glass bottles 1 L or 2 L on OMNIS Reagent Organizer.



6.02050.110

Bottle fixation RO, universal

Bottle fixation for glass bottles 1 L or 2 L on OMNIS Reagent Organizer.



6.02765.000

Waste collector, 4 x UNF 10 /32

Waste collector for up to four waste capillaries.



6.05332.000

OMNIS IC equip.: Container 10 L, LC

Accessory set for automatic monitoring of 10 L container.

Prerequisites:

OMNIS IC with ROLC



6.05801.040

Volume sensor, ROLC kit

Upgrade kit to extend the OMNIS Reagent Organizer Liquid Control with an additional volume sensor.

Prerequisites:

– 2.1056.1200 OMNIS Reagent Organizer Liquid Control

Info: This upgrade must be installed by an authorized Metrohm service technician.



6.05801.090

Magnetic stirrer, ROLC kit



Upgrade kit to extend the OMNIS Reagent Organizer Liquid Control with a magnetic stirrer.

Prerequisites:

– 2.1056.1200 OMNIS Reagent Organizer Liquid Control

Info: This upgrade must be installed by an authorized Metrohm service technician.
