

Pump Unit 772 "Rinse"

2.772.0030

Peristaltic pump for rinsing sample beakers and electrodes. The "Rinse" version contains extensive accessories for rinsing. The pump can be actuated through a Relay Box 731. It is suitable for the conveyance of solvents.

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 2.772.0030

Order Number	Description	Quantity
6.1621.000	Container 10 L As rinsing or waste container in automated systems.	1 PCS
		
6.1805.060	FEP tubing / M6 / 60 cm With light and kink protection	3 PCS
		
6.1812.000	PTFE tubing 4 mm / 6 mm / 4 m	1 PCS
		

6.1820.050

M6 / Olive adapter

2 PCS



For 6.1826.100 pump tubing, to peristaltic pumps.

6.1826.100

Pump tubing 6.4 mm / 1 m

1 PCS



Replaces the pump tubing used with the 772 Pump Unit and the 843 Pump Station (peristaltic).

6.1828.000

PVDF connection nipple

1 PCS



For 6.1621.000 canister.

6.2740.020

Spray nozzle

3 PCS



For the fine-spraying of the rinsing solution.

Optional accessories

6.1826.160

Viton pump tubing

For peristaltic pumps.

