



786 Swing Head, right-swinging, reinforced

2.786.0240

Strengthened Swing Head for usage with Sample Processors.

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 2.786.0240

Order Number	Description	Quantity
1.786.0010	786 Swing Head with transfer head left Additional equipment for Robotic Sample Processors XL for transferring / pipetting from small sample vessels into larger titration vessels on the rack or into an external titration cell; this can be located to the left of the work station.	1 PCS
6.2058.090	Reinforcement for SwingHead R Reinforcement for 786 Swing Head (series 31 and higher) for usage of swing arms with stripping function	1 PCS
6.2621.120	Hexagon key 1.5 mm	1 PCS

Optional accessories

6.1462.240

**Transfer head for 786 Swing
Head, bent and right
swinging**



Transfer head, bent and right swinging, for 786 Swing Head. In combination with the 6.1808.220 Adapter the transfer head can be used on multi-row racks to pick up tools with luer connection.