




801 Stirrer

2.801.0010

Magnetic stirrer without stand for supplementing the Titrino plus, Dosimat plus, Titrandos, Sample Processors, 805 Dosimats and 780/781 pH Meters. With permanently attached cable for MSB (Metrohm Serial Bus).

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 2.801.0010

Qt.	Order no.	Description
1 PCS	1.801.0010	801 Stirrer
Magnetic stirrer without stand for supplementing the Titrino plus, Dosimat plus, Titrandos, Sample Processors, 805 Dosimats and 780/781 pH Meters. With permanently attached cable for MSB (Metrohm Serial Bus).		



1 PCS

6.1903.020

Stirring bar / 16 mm

Stirring Bar with magnetic core, PTFE coated, length 16 mm



1 PCS

6.1903.030

Stirring bar / 25 mm

Stirring bar with magnetic core, PTFE covering, length 25 mm.



Optional accessories

Order no.	Description	
6.2001.090	Base plate on both sides Double stand for attaching 2 x 801, 803, 804 to a Titrande	
6.2151.010	Cable MSB plug / socket Extension cable for MSB connections and for Touch Control - Titrande, Mini-DIN socket - mini-DIN plug	
6.2070.000	Base enlargement to stirrers 80X Enables to stir 1 and 2 L bottles on the 80X stirrers.	