



Comb. iAg ring electrode

6.0450.300

Intelligent, combined silver ring electrode with integrated memory chip for storing of sensor data.
This electrode is suitable for precipitation titrations when the pH value varies (titrant silver nitrate) of, e.g.:

- Chloride, Bromide, Iodide
- Sulfides
- Hydrogen sulfide
- Mercaptans
- Cyanides

$c(\text{KNO}_3) = 1 \text{ mol/L}$ is used as reference electrolyte and for storage.

iTrodes can be used on Titrando, Ti-Touch or 913/914 meters.

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 6.0450.300

Order Number	Description	Quantity
6.1236.050	Sleeve with SGJ 14/12 mm Sleeve with SGJ 14/12 mm, polyethylene.	1 PCS
		
6.2008.040	Storage vessel Together with 6.2008.050 storage vessel holder. it provides a support for the electrode on 807 Dosing Units.	1 PCS
		

Optional accessories

6.2310.010

Electrolyte KNO₃-1mol/L 250 mL

Electrolyte solution KNO₃ 1 M (reference electrolyte for combined silver electrodes)



2.854.0010

854 iConnect

The 854 iConnect – Electrode cable and measuring amplifier for intelligent "iTrodes" electrodes (cable length 150 cm).

