




Accessory kit: OMNIS IC L Module

6.05000.040

Accessory set to complete the installation of the OMNIS IC L module.

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 6.05000.040

Order Number	Description	Quantity
6.02744.010	Flangeless fitting M6, 1/8", 2 pcs Flangeless fitting for 1/8" tubing. Quantity: 2 pcs	2 PCS
		
6.02837.040	Filter tube Filter tube with Luer adapter, protects eluents and reagents from particles.	2 PCS
		
6.01806.000	PEEK capillary, 2000/0.18 mm	1 PCS
		

6.01805.500

FEP tubing 1/8" RO – DEG, 1,200 mm

2 PCS

Pre-cut FEP tubing with identification label for OMNIS IC



6.01831.100

PEEK capillary DEG – HPP, 250/0.75 mm

2 PCS

Pre-cut PEEK capillary with identification label for OMNIS IC



6.01831.110

PEEK capillary HPP – INJ, 160/0.5 mm

2 PCS

Pre-cut PEEK capillary with identification label for OMNIS IC




6.1825.210

PEEK sample loop 20 µL

2 PCS

For injection valve, with 2 PEEK pressure screws



6.01831.140

PEEK capillary INJ – COL, 180/0.18 mm

1 PCS

Pre-cut PEEK capillary with identification label for OMNIS IC



6.2744.070

Pressure screw short

2 PCS

Short version. With UNF 10/32 connection. 5 pieces. For the connection of PEEK capillaries



6.2744.320

PEEK coupling M6 / UNF 10/32

2 PCS

Connecting transfer tubing (M6) with capillaries (UNF 10/32).



6.01831.180

PEEK capillary INJ – COL, 280/0.18 mm

1 PCS

Pre-cut PEEK capillary with identification label for OMNIS IC



6.05801.310

Drain kit, L

1 PCS



Drain kit to replace the drain tubing and connections in the OMNIS IC L.

Prerequisites:

2105200x0 OMNIS IC L