



NH₃-selective gas membrane electrode (high conc.)

6.0506.150

Combined ammonia-selective gas membrane electrode for higher concentrations (compared to 6.0506.100), e.g.:

- in waste water
- in samples containing oil
- for long-term measurements

An individual contaminated membrane can easily be cost-effectively replaced.

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 6.0506.150

Order Number	Description	Quantity
6.1236.050	Sleeve with SGJ 14/12 mm Sleeve with SGJ 14/12 mm, polyethylene.	1 Stk
		
6.2008.040	Storage vessel Together with 6.2008.050 storage vessel holder. it provides a support for the electrode on 807 Dosing Units.	1 Stk
		

Optional accessories

6.2316.030

**Inner filling solution to
6.0506.100 electrode**

Measuring electrolyte for NH_3 -selective electrode (50 mL)



6.2104.020

Electrode cable / 1 m / F

For connecting electrodes with Metrohm plug-in head G to
Metrohm instruments (socket F).

