



KFC equipment for OMNIS Coulometer with generator electrode with diaphragm

6.05406.000

Complete accessory set for coulometric Karl Fischer titration.

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 6.05406.000

Order Number	Description	Quantity
6.00348.100	Generator electrode with diaphragm  <p>Generator electrode with diaphragm for coulometric Karl Fischer titrations, SGJ 29/22. This electrode has a cable connection angled upwards, which simplifies cable guiding and operation when used with the OMNIS Sample Robot Oven.</p> <p>This generator electrode is well-suited for:</p> <ul style="list-style-type: none"> – Samples with very low absolute water content 50 ppm – Coulometric water content determination in oil samples – Applications where ketone reagents are used – Applications where low-conducting reagents are used (addition of chloroform, xylene >10%) 	1 PCS
6.02047.000	OMNIS titration cell holder  <p>Titration cell holder for fastening a coulometric titration cell to an OMNIS Coulometer or OMNIS Coulometer Module</p>	1 PCS

6.02104.040

Electrode cable plug-in head G (pol.) / plug P, 0.55 m 1 PCS

Cable for connecting an analog polarizable electrode (plug-in head G) to an analog measuring module from OMNIS.



6.0341.100

Double Pt-wire electrode for coulometry 1 PCS

Indicator electrode used for coulometric Karl Fischer titration.



6.1403.030

KF adsorber tube for coulometer cell 1 PCS

For KF coulometers.



6.1437.000

SGJ stopper SGJ 14 1 PCS


6.1448.020

Septum 16 mm, 5 pieces

2 PCS

Set of 5 pcs.



6.1464.320

KF titration vessel / 80-250 mL / coulometric

1 PCS



6.1903.030

Stirring bar / 25 mm

1 PCS

Stirring bar with magnetic core, PTFE covering, length 25 mm.



6.2104.620

Electrode cable / 0.65 m / H

1 PCS

With plug H, for the generator electrode - coulometer connection



6.2701.040



GL 18 screw cap with hole

1 PCS

Made of PBT. With a hole for the 6.1448.020 septum. For the 6.5405.000 cell without diaphragm and the 6.1464.32 X and 6.1465.320 titration vessels.

6.2713.000



PTFE sleeve SGJ14

2 PCS

For fat-free ground-joint connection with the standard ground-joint SGJ 14.

6.2713.010



PTFE sleeve SGJ29

1 PCS

Sleeve for fat-free ground-joint connection with the standard ground-joint SGJ 29.

6.2713.020



PTFE sleeve SGJ19

1 PCS

For fat-free ground-joint connection with the standard ground-joint SGJ 19.

Optional accessories

6.2811.000

Molecular sieve



Molecular sieve. Bottle containing 250 g. Pore size: 0.3 nm. Without moisture indicator. For Rancimats and Karl Fischer instruments.

6.2811.010

Molecular sieve (US)



Molecular sieve. Bottle containing 250 g. Pore size: 0.3 nm. Without moisture indicator. For Rancimats and Karl Fischer instruments.

This product is intended for all users who are US residents. The material is rated hazardous in the US in accordance with OSHA Hazard Communication Standard (29 CFR 1910.1200).