



Metrosep C 6 - 150/4.0

6.1051.420

The high-capacity C 6 material makes the Metrosep C 6 - 150/4.0 separation column the optimum solution for separating standard cations with high differences in concentration in conjunction with reasonable retention times. Drinking water with low ammonium contents can be determined with this column.

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 6.1051.420

Order Number	Description	Quantity
6.2744.060	Threaded stopper For UNF 10/32. Stopper for IC, e.g. for the sealing of columns.	2 PCS



Optional accessories

6.1051.500

Metrosep C 6 Guard/4.0



The Metrosep C 6 Guard/4.0 contains the C 6 column material and is used to protect against particles and contamination. This considerably lengthens the service life of the analytical separation column. The Metrosep C 6 Guard/4.0 works according to the "On Column Guard System" and is screwed directly onto the respective separation column with virtually no dead volume.

6.1051.510

Metrosep C 6 S-Guard/4.0



The Metrosep C 6 S-Guard/4.0 contains the C 6 column material and is used to protect it against particles and contamination. This considerably lengthens the service life of the analytical separation column.