



Fastening plate for deflector and holder for collection container

6.2058.050

The fastening plate is mounted at the side of the Sample Processor tower and allows the mounting of a deflector for the automatic pulling off the used tools.

Below, the accessories are grouped into Scope of delivery and Optional accessories. Please keep this printout at hand for ordering replacement material. These lists may be subject to change.

Optional accessories

Order no.	Description
-----------	-------------

6.2057.150

Holder for collection container

Holder for mounting collection container



6.2058.060

Deflector for adapter 61808220

The deflector is mounted with the fastening plate (6.2058.050) at the side of the Sample Processor tower and allows the automatic pulling off the used luer tools into the collection container (6.1625.010)



6.2058.070

Deflector for adapter 61808250

The deflector is mounted with the fastening plate 6.2058.050 at the side of the Sample Processor tower and allows the automatic pulling off the used pipetting tips into the collection container 6.1625.010.

