




Sample needle made of steel, coated

6.2624.200

Alternative to the sample needle made of zirconium oxide. Particularly suitable for septum vessels. Metal-free thanks to special coating. Requires the 1/16 in. needle holder (6.2833.010) for installation on the Sample Processor.

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Optional accessories

6.2629.000

Screw nut to needle holder IC

On Sample Processors



6.2833.010

Needle/tube holder 1/16 in.

For IC Sample Processors. For 1/16 in. tubes

