



VA electrode equipment with Multi-Mode Electrode pro for Professional VA instruments

6.5339.030

Complete electrode set for polarographic and voltammetric determinations. Includes Multi-Mode Electrode pro, reference electrode, platinum auxiliary electrode, measuring vessel, stirrer, electrolyte solution and additional accessories for setting up and operating the Multi-Mode Electrode.

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 6.5339.030

Order Number	Description	Quantity
6.0343.100	Separate Pt Rod Electrode Platinum auxiliary electrode for voltammetry and CVS.	1 PCS
		
6.0728.120	Ag/AgCl reference electrode Inner reference system with electrolyte $c(\text{KCl}) = 3 \text{ mol/L}$. For use with electrolyte vessel made of plastic 6.1245.010.	1 PCS
		
6.1204.500	Stirrer for Professional VA instruments Stirrer for Professional VA instruments, complete, metal-free.	1 PCS
		

6.1226.030


Mercury drop capillary (non-silanized) for MME

1 PCS

Non-silanized mercury drop capillary (10 pieces)

6.1226.050


Mercury drop capillary (silanized) for MME

1 PCS

Silanized mercury drop capillary (10 pieces)

6.1245.010


Electrolyte vessel made of plastic

1 PCS

Electrolyte vessel made of plastic with ceramic diaphragm for double junction reference electrode.

6.1246.120


Multi-Mode Electrode pro

1 PCS

Mercury electrode for voltammetry. Can be operated as DME, SMDE or HMDE.

6.1247.020

Needles for MME

1 PCS



Needles for Multi Mode Electrode (3 pieces).

6.1415.210

Titration vessel / 10-90 mL

1 PCS



6.2308.020

Electrolyte 3 mol/L KCl (250 mL)

1 PCS

Electrolyte solution $c(\text{KCl}) = 3 \text{ mol/L}$, (for Ag/AgCl reference systems)


6.2406.000

Mercury drop catcher

1 PCS



Silver wire (to capture mercury drops by amalgamation)

6.2711.030

Tray for filling the MME

1 PCS

For 884 Professional VA, 797 VA Computrace, 757 VA Computrace and 647, 663, 694, 747 VA Stands



6.2816.020

Syringe 10 mL with Luer connection

1 PCS

For various applications in IC and VA



6.2816.030

Needle with Luer connector

1 PCS

Needle for filling in mercury into the Multi-Mode Electrode and for sample addition in Karl Fischer titration.

