




KFT equipment for 810/814 Sample Processor

6.5610.010

Karl Fischer equipment for 810 or 814 Sample Processor.

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 6.5610.010

Order Number	Description	Quantity
6.0340.000	Double Pt-wire electrode for sample changer (fixed cable) Indicator electrode (diameter = 5.3 mm, fixed cable 2 m) used for automated Karl Fischer titration.	1 PCS
		
6.1432.210	Sample beaker, clear glass, 75 mL, 1 piece Sample beaker, clear glass, 75 mL, suitable for sample racks: 6.2041.340, 6.2041.350, 6.2041.800 and 6.2041.900, 1 piece	24 PCS
		
6.1458.030	KF titration head / 4 x M10 connectors With 4 M10 screw nipples, including stoppers.	1 PCS
		

6.1805.060

FEP tubing / M6 / 60 cm

1 PCS

With light and kink protection



6.1808.000

Connection bushing M6

1 PCS



6.1820.020

M10/M6 adapter

1 PCS

Connector with M10 thread for 6.1812.XXX tubing and M6 thread for 6.1805.XXX tubing.



6.1821.000

Aspiration tube

1 PCS

For KF Sampler



6.1903.010

Stirring bar / 12 mm

24 PCS

Stirring Bar with magnetic core, PTFE coated, length 12 mm



6.2037.040

Foil holder

24 PCS

Foil holder for 75 mL sample beaker for Karl Fischer for KF Sample Changer



6.2041.340

Sample rack, 24 x 75 mL, (PVC)

1 PCS

Sample rack for 810 Sample Processor, 814 USB Sample Processor, 855 Robotic Titrosampler "Light", suitable for 24 sample beakers. The following sample beakers can be used: 6.1432.210, 6.1459.400.
Plastic: Polyvinyl chloride (PVC)



6.2820.000

Al foil

1 PCS

1,000 rondels with a diameter of 80 mm. Foil thickness: 0.010 mm. Together with 6.2037.040 foil holder for 6.1432.210 sample beaker and 6.2037.080 foil holder for 6.1459.300 sample beaker.



Optional accessories

6.2811.000

Molecular sieve



Molecular sieve. Bottle containing 250 g. Pore size: 0.3 nm. Without moisture indicator. For Rancimats and Karl Fischer instruments.

6.2811.010

Molecular sieve (US)



Molecular sieve. Bottle containing 250 g. Pore size: 0.3 nm. Without moisture indicator. For Rancimats and Karl Fischer instruments.

This product is intended for all users who are US residents. The material is rated hazardous in the US in accordance with OSHA Hazard Communication Standard (29 CFR 1910.1200).