




Consumable Kit CVS

6.5703.020

Assembly of important expendable items for CVS applications with Professional CVS instruments.

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 6.5703.020

Order Number	Description	Quantity
6.0728.130	LL-Ag/AgCl reference electrode for CVS Inner reference system with electrolyte $c(\text{KCl}) = 3 \text{ mol/L}$. For use with electrolyte vessel made of plastic 6.1245.010.	1 PCS
		
6.1204.610	2 mm platinum electrode tip for CVS Metal electrode tip made of platinum (Pt), diameter of electrode disk $2 \text{ mm} \pm 0.05 \text{ mm}$, polished surface, shaft made of glass, M3 thread connection. Concentricity error: 0.3 mm. Temperature range: 0 - 50 °C.	1 PCS
		
6.1244.050	Drive belts for rotating disk electrodes For Professional VA/CVS instruments (3 pcs.).	1 PCS
		

6.1824.000

Four-way micro dosing tip

1 PCS

With 4 x M6 thread. For VA measuring stands.

