





Liquid Autosampler CIC (TEI), 30 positions

2.0136.0610


The Liquid Autosampler CIC (TEI) is used for automated analysis of liquid samples using Combustion IC. It can be used in combination with the Boat Introduction Modules (6.07311.010 or 6.07311.020) and the Liquid Introduction Module (6.07311.040).

Phạm vi giao hàng 2.0136.0610

Order Number	Description	Quantity
1.0136.0610	Liquid Autosampler CIC (TEI), 30 positions	1 PCS



6.7303.110	Snap cap with red septum, 100 pieces	1 PCS
------------	---	-------



Snap cap with septum for 11 mm vials. Provides reliable sealing and broad chemical compatibility for routine laboratory analysis.

Quantity: 100 pieces

6.7303.100	Vials glass, 1.8 ml, U-shape, 100 pieces	1 PCS
------------	---	-------



Maximize sample recovery with this universal 1.8 mL snap ring vial designed for low-volume analytical applications.

Quantity: 100 pieces

6.07313.040

Wash/waste container set for Liquid Autosampler CIC (TEI) 1 PCS


Wash/waste container set for the Liquid Autosampler CIC (TEI).

6.07311.290

Septum for Introduction Module (TEI) 1 PCS


Septum for liquid injection port of the Introduction Modules (TEI) 6.07311.010, 6.07311.020, 6.07311.030, 6.07311.040.

6.07313.020

Screw incl. ferrule for Liquid Autosampler CIC (TEI) 1 PCS


Screw including collapsible ferrule for 1/16" tubing to connect the transfer line to the needle on the Liquid Autosampler CIC (TEI).

6.07313.110

Replacement needle for CIC (TEI) 1 PCS


Replacement needle (TEI) for dosing syringe for Liquid Autosampler CIC (2.0136.0610) or manual liquid sample injection.

6.07313.050

Prime Vial Set for Liquid Autosampler CIC (TEI)

1 PCS

Set of Prime Vials for the preparation of the Liquid Autosampler CIC (TEI).



Phụ kiện tùy chọn

6.07313.020

Screw incl. ferrule for Liquid Autosampler CIC (TEI)

Screw including collapsible ferrule for 1/16" tubing to connect the transfer line to the needle on the Liquid Autosampler CIC (TEI).



6.07313.100

Syringe 100 µl manual injection (TEI)

Dosing syringe for manual injection of liquid samples using one of the three Introduction Modules (6.07311.010, 6.07311.020, or 6.07311.030).



6.07313.110

Replacement needle for CIC (TEI)

Replacement needle (TEI) for dosing syringe for Liquid Autosampler CIC (2.0136.0610) or manual liquid sample injection.



6.7302.120

**Quartz wool for 6.7302.110,
10 g**

Quartz wool for fixing the sample in place in the sample vial
made of quartz (6.7302.110).

