



OMNIS Dosing Module without stirrer

2.1003.0010

Dosing module for connection to an OMNIS Titrator for extending the system to include an additional buret for titration/dosing. Can be supplemented with one magnetic stirrer or rod stirrer for use as separate titration stand. Freely selectable cylinder unit with 5, 10, 20 or 50 mL.

Phạm vi giao hàng 2.1003.0010

Order Number	Description	Quantity
6.01600.010	OMNIS Liquid Adapter Adapter for connecting reagent bottles with 3S technology to an OMNIS Titrator, Titration Module or Dosing Module.	1 PCS
6.01601.000	Bottle cap, multi-use Reusable bottle cap with 3S technology for connecting reagent bottles to an OMNIS Titrator, Titration Module or Dosing Module.	1 PCS
6.1608.020	Glass bottle / 1000 mL / GL 45 Bottle for auxiliary solutions.	1 PCS



6.1546.030

**Piston tongs**

1 PCS

For pistons for 807 Dosing Units or cylinder units OMNIS.

6.02004.010

**Bottle holder**

1 PCS

Bottle holder for OMNIS Titrator, Titration Module and Dosing Module.

Phụ kiện tùy chọn

6.02102.010
Cable MDL PL/SO 0.5 m

Connecting cable for OMNIS instruments with MDL connector
(0.5 m)



6.02102.020
Cable MDL PL/SO 1 m

Connecting cable for OMNIS instruments with MDL connector (1
m)



6.02102.030
Cable MDL PL/SO 2 m

Connecting cable for OMNIS instruments with MDL connector (2
m)



6.02102.040

Cable MDL PL/SO 5 m

Connecting cable for OMNIS instruments with MDL connector (5 m)



6.02102.050

Cable MDL PL/SO 10 m

Connecting cable for OMNIS instruments with MDL connector (10 m)



6.03001.120

Cylinder unit OMNIS 2 mL

Intelligent 2 mL cylinder unit for an OMNIS Titrator, Titration Module or Dosing Module. Includes dosing tubing and antidiffusion tip.



6.03001.150

Cylinder unit OMNIS 5 mL

Intelligent 5 mL cylinder unit for an OMNIS Titrator, Titration Module or Dosing Module. Includes dosing tubing and antidiffusion tip.



6.03001.210

Cylinder unit OMNIS 10 mL

Intelligent 10 mL cylinder unit for an OMNIS Titrator, Titration Module or Dosing Module. Includes dosing tubing and antidiffusion tip.



6.03001.220

Cylinder unit OMNIS 20 mL

Intelligent 20 mL cylinder unit for an OMNIS Titrator, Titration Module or Dosing Module. Includes dosing tubing and antidiffusion tip.



6.03001.250

Cylinder unit OMNIS 50 mL

Intelligent 50 mL cylinder unit for an OMNIS Titrator, Titration Module or Dosing Module. Includes dosing tubing and antidiffusion tip.



6.01503.120

**Cylinder unit OMNIS 2 mL
without accessories**

Intelligent cylinder unit OMNIS 2 mL with four M6 ports for OMNIS Titrator, Titration Module, Dosing Module or OMNIS Dosing Drive Unit.



6.01503.150

**Cylinder unit OMNIS 5 mL
without accessories**

Intelligent cylinder unit OMNIS 5 mL with four M6 ports for OMNIS Titrator, Titration Module, Dosing Module or OMNIS Dosing Drive Unit.



6.01503.210

**Cylinder unit OMNIS 10 mL
without accessories**

Intelligent cylinder unit OMNIS 10 mL with four M6 ports for OMNIS Titrator, Titration Module, Dosing Module or OMNIS Dosing Drive Unit.



6.01503.220

**Cylinder unit OMNIS 20 mL
without accessories**

Intelligent cylinder unit 20 mL for one OMNIS Titrator, Titration Module or Dosing Module. Without dosing tubing and antidiffusion tip.



6.01503.250

**Cylinder unit OMNIS 50 mL
without accessories**

Intelligent cylinder unit OMNIS 50 mL with four M6 ports for OMNIS Titrator, Titration Module, Dosing Module or OMNIS Dosing Drive Unit.



6.02106.000

**Liquid Adapter extension
cable, 0.55 m**

Extension cable (0.55 m in length) for connecting an OMNIS Liquid Adapter to an OMNIS Titrator, Titration Module or Dosing Module.



6.02106.010

**Liquid Adapter extension
cable, 2 m**

Extension cable (2 m in length) for connecting an OMNIS Liquid Adapter to an OMNIS Titrator, Titration Module or Dosing Module.



6.05609.000

**KFT equipment for OMNIS
Titrator**

Complete accessory set for volumetric Karl Fischer titration.



6.05609.010

**KFT equipment for OMNIS
Titrator with Polytron**

Complete accessory set for volumetric Karl Fischer titration with Polytron.

