



Eco KF Titrator

2.1027.0010

The new Eco KF Titrator with integrated magnetic stirrer und touch-sensitive display is ideal for routine water content determination. Predefined methods enable a smooth and easy instrument start-up. In addition to GLP-compliant printouts on paper or as PDFs, the Eco KF Titrator also offers the option of connecting balances or of sending the determination data to a PC via PC / LIMS reports.

Phạm vi giao hàng 2.1027.0010

Qt.	Order no.	Description
1 PCS	1.1027.0010	Eco KF Titrator

The new Eco KF Titrator with integrated magnetic stirrer und touch-sensitive display is ideal for routine water content determination. Predefined methods enable a smooth and easy instrument start-up. In addition to GLP-compliant printouts on paper or as PDFs, the Eco KF Titrator also offers the option of connecting balances or of sending the determination data to a PC via PC / LIMS reports.



1 PCS

6.1446.080

Threaded stopper / M8



1 PCS

6.1546.040

Piston tool for Eco Dosimat/Eco Titrator

Tool for removing and mounting the piston.



1 PCS

6.1602.105

Bottle attachment / GL 45 / green

For direct dosing from reagent bottles with GL 45 thread.



1 PCS

6.1619.010

Adsorber complete, for Exchange Unit



2 PCS 6.1805.090 FEP tubing / M6 / 31 cm
With light and kink protection



1 PCS 6.1805.100 FEP tubing / M6 / 40 cm
With light and kink protection



1 PCS 6.1819.020 Tube
For M6 threaded holes. For 6.1602.1XX bottle attachments



1 PCS 6.2013.010 Clamping ring
For support rods with a diameter of 10 mm.



1 PCS

6.2016.070

Support rod / 400 mm



1 PCS

6.2023.020

Clip for SGJ 14/15

Clip for SGJ 14/16



1 PCS

6.2164.010

Power supply unit 100 - 240 V/24 V DC

Power supply unit for 915 KF Ti-Touch, 916 Ti-Touch, 917 Coulometer and all Eco Titrators / Eco Dosimat.



1 PCS

6.02712.000

Protective cover for Eco Dosimat / Eco Titrator

Protective cover for Eco Dosimat / Eco Titrator display.



1 PCS

6.6081.020

USB flash drive for Eco KF Titrator

Storage medium for Eco KF Titrator



Phụ kiện tùy chọn

Order no.	Description	
2.1029.0010	Solvent Pump	
	Solvent Pump is used with an Eco Titrator / Eco KF Titrator / Eco Coulometer for reliably adding solvent and aspirating off titration vessel content without contact with the reagents. Device-independent use is also possible with optional 6.2164.010 power supply unit.	
6.02709.040	Stopper	
	Stopper for titration vessel lids 6014050X0.	
6.02711.000	Spoon for paste	
	Suitable for use in Karl Fischer water content determinations of pasty substances such as honey, cosmetics, greases or various dairy products. The spoon for paste facilitates the addition of viscous samples. It is handled like any other spatula and enables convenient, rapid, and loss-free addition of pastes into the Karl Fischer titration vessel. We furthermore recommend the use of a homogenizer for dissolving viscous samples quickly and completely. Compatible with 6014050X0 titration vessel lids.	
6.05610.050	KFT equipment for Eco KF Titrator	
	Complete accessory set for volumetric Karl Fischer titration	

6.05610.060 KFT equipment for Eco KF Titrator with Polytron

Complete accessory set for volumetric Karl Fischer titration with Polytron



6.03002.150 Cylinder unit Eco 5 mL

Cylinder unit 5 mL for Eco Dosimat / Eco Titrator. Without dosing tubing and antidiffusion tip.



6.03002.210 Cylinder unit Eco 10 mL

Cylinder unit 10 mL for Eco Dosimat / Eco Titrator. Without dosing tubing and antidiffusion tip.



6.03002.220 Cylinder unit Eco 20 mL

Cylinder unit 20 mL for Eco Dosimat / Eco Titrator. Without dosing tubing and antidiffusion tip.



6.03002.250 Cylinder unit Eco 50 mL

Cylinder unit 50 mL for Eco Dosimat / Eco Titrator. Without dosing tubing and antidiffusion tip.



2.142.0100 Custom Q3X thermal printer

Compact printer with USB interface for

- 900 Touch Control
- 915 KF Ti-Touch
- 916 Ti-Touch
- 917 Coulometer
- 877 / 848 Titrino plus
- 865 / 876 Dosimat plus
- 91X Meter (cable 6.2151.140)
- Eco Dosimat / Titrator
- 862 Compact Titrosampler
- 870 KF Titrino plus
- 899 Coulometer



Paper width 60 mm (40 characters). Including 6.2151.120 USB cable.

6.9001.050 Thermal paper for Custom Neos-N / Q3X

Thermal paper 60 mm / 25 m for Custom Neos-N / Q3X.
New: 25-year legibility with correct storage.

1 roll



Compact impact printer with USB and RS232 interface for

- 900 Touch Control (+)
- 915 KF Ti-Touch (+)
- 916 Ti-Touch (+)
- 917 Coulometer (+)
- 91X Meter (+)
- Eco Dosimat/Titrator (+)
- 877/848 Titrino plus (additional 6.2151.100 required)
- 865/876 Dosimat plus (additional 6.2151.100 required)
- 862 Compact Titrosampler (additional 6.2151.100 required)
- 870 KF Titrino plus (additional 6.2151.100 required)
- 899 Coulometer (additional 6.2151.100 required)



(+) Cable included in scope of delivery

6.01406.250

50 – 150 mL titration vessel

50 – 150 mL titration vessel with screw top for use with 6014050X0 titration vessel lids.



6.2811.000

Molecular sieve

Molecular sieve. Bottle containing 250 g. Pore size: 0.3 nm. Without moisture indicator. For Rancimats and Karl Fischer instruments.



6.2811.010

Molecular sieve (US)

Molecular sieve. Bottle containing 250 g. Pore size: 0.3 nm. Without moisture indicator. For Rancimats and Karl Fischer instruments.

This product is intended for all users who are US residents. The material is rated hazardous in the US in accordance with OSHA Hazard Communication Standard (29 CFR 1910.1200).

