




786 Swing Head with transfer head right

2.786.0020

Additional equipment for Robotic Sample Processors XL for transferring / pipetting from small sample vessels into larger titration vessels on the rack or into an external titration cell; this can be located to the right of the work station.

Phạm vi giao hàng 2.786.0020

Order Number	Description	Quantity
1.786.0010	786 Swing Head with transfer head left Additional equipment for Robotic Sample Processors XL for transferring / pipetting from small sample vessels into larger titration vessels on the rack or into an external titration cell; this can be located to the left of the work station.	1 PCS
6.1562.100	Pipetting tubing / 10 mL Pipetting tubing for fully automatic transfer of sample solutions.	1 PCS
6.2621.120	Hexagon key 1.5 mm	1 PCS



6.1462.040

Robotic arm with transfer head for 786 Swing Head, right swinging 1 PCS

Transfer robotic arm for liquid handling purposes



6.1823.010

Guide shaft to pipetting tubing 1 PCS

For assembly in the transfer robotic arm on the Robotic Sample Processors

