



Combined dCa ISE

6.00502.300

Digital, combined calcium-selective electrode for OMNIS.

This ISE is suitable for:

- Ion measurements of Ca^{2+} in aqueous solutions
- Complexometric (back) titrations (e.g., determination of water hardness)

Thanks to a robust/break-proof plastic shaft made of propylene and an impact protection for the polymer membrane, this sensor is mechanically very resistant.

The reference electrolyte used is $c(\text{KCl}) = 3 \text{ mol/L}$.

dTrodes can be used on OMNIS Titrators.

Phạm vi giao hàng 6.00502.300

Order Number	Description	Quantity
6.1236.050	Sleeve with SGJ 14/12 mm	1 Stk



Sleeve with SGJ 14/12 mm, polyethylene.

6.2008.040

Storage vessel

1 Stk



Together with 6.2008.050 storage vessel holder. it provides a support for the electrode on 807 Dosing Units.

Phụ kiện tùy chọn

6.2308.020

Electrolyte 3 mol/L KCl (250 mL)

Electrolyte solution $c(\text{KCl}) = 3 \text{ mol/L}$, (for Ag/AgCl reference systems)



6.02104.300

Electrode cable plug-in head Q / plug P, 0.55 m

Cable for connecting a dTrobe (plug-in head Q) to the digital measuring module from OMNIS.



6.02104.310

Electrode cable plug-in head Q / plug P, 1.5 m

Cable for connecting a dTrobe (plug-in head Q) to the digital measuring module from OMNIS.

